



Changes in the Detection of Clinically Relevant White Matter Lesions in MS Using FDA-Cleared Automated Software

Kelly Leyden MRes¹, Anisha Keshavan PhD¹, Victor Gehman PhD¹, Michael Iv MD^{1,2}

¹Octave Bioscience, Menlo Park, CA, United States of America

²Stanford University, CA, United States of America

DISCLOSURE



Kelly Leyden, MRes

Product Lead - MRI Insights
Employee of **Octave Bioscience**

Anisha Keshavan, PhD

Senior Data Scientist
Employee of **Octave Bioscience**



Michael Iv, MD

Neuroradiologist
Octave Bioscience
Clinical Associate Professor of Radiology
Stanford University

Contribution to this study was as a paid consultant, and was not part of his Stanford University duties or responsibilities

Victor Gehman, PhD


Principal Data Scientist
Employee of **Octave Bioscience**



BACKGROUND

Improving lesion burden detection in patients with multiple sclerosis (MS)

- Outside of gadolinium enhancing T1 lesions, disease activity in MS is often characterized by new, enlarging, and shrinking T2 lesions on MRI.
- Use of institution-specific semi-automated software can improve white matter lesion detection, which can influence neurologists' clinical decision-making.¹
- Commercially available FDA-cleared automated software can potentially provide similar benefits but with the added advantage of being widely available in the United States.




¹ Van Heerden J et al. *AJNR*, 2015

OBJECTIVE


Improved white matter lesion detection with FDA-cleared automated software?

- To compare and characterize the detection of **clinically relevant** white matter lesions on MRI without and with the use of a commercially available FDA-cleared automated software package.


New T2 Lesion



Enlarging T2 Lesion




Shrinking T2 Lesion




APPROACH

Image acquisition and post-processing

- Retrospective study w/ subject selection from anonymized dataset between 2013-2019
- 2 MRIs/patient: 3D T1 and 3D T2 FLAIR
- 3D T1 and 3D T2 FLAIR images processed using FDA-cleared software (LesionQuant 3.0.1) to produce:
 - **Co-registered, color-coded map highlighting lesion changes between MRI comparisons**
 - **Hot (red/yellow): new or enlarging lesion**
 - **Cold (blue): shrinking lesion**



Comparison of radiology findings



- Using only visual interpretation, a board-certified neuroradiologist reported the number of new, enlarging, and shrinking T2 lesions.
- To avoid recall bias, the second interpretation with the aid of the software package (including processed data and the color-coded head map) was performed after a one-month waiting period.
- Lesion differences were compared with paired t-tests.
- Agreement was assessed with Bland-Altman analysis.

RESULTS

Detection of MS lesions without and with software use

- 52 MS patients who had 2 brain MRIs (mean MRI time interval: 12 months)
- In **23 patients (44%)**, new or enlarging lesions were found with the software that were missed or misidentified on visual inspection alone, yielding on average **2 more new or enlarging lesions per patient**.
- **Good agreement** between both methods, but the use of software resulted in **detection of more lesions in patients with high lesion burden** (mean difference: 0.7; 95% limits of agreement: -2.9, 4.2)


Comparison of Lesion Counts

	W/O software	W/ software	% difference
New	67	82	20.1% P=0.11
Enlarging	13	32	84.4% P=0.02*
Shrinking	14	15	6.9% P=0.57


RESULTS

Examples:
Lesions
missed with
visual
inspection
only


***New
T2 Lesion***



***Enlarging
T2 Lesion***



***Shrinking
T2 Lesion***



DISCUSSION

Software detects more clinically relevant MS lesions

- Use of FDA-cleared automated software improved the detection of MS WMLs, especially of enlarging T2 lesions in patients with **high lesion burden**.
- In nearly half of patients, the number of T2 lesion differences detected with the software was different from the routine evaluation, suggesting an **opportunity to improve reporting of disease activity and impact clinical outcomes**.

NEXT STEPS

Assess impact on clinical decision-making

- Evaluate additional metrics such as T1 WM hypointensities and brain volumes
- Assess MS Neurologists' satisfaction and assess potential change in treatment decision-making based on improved neuroimaging findings

Questions or comments?

Corresponding Author: kleyden@octavebio.com

Visit additional posters in collaboration with Octave Bioscience

P010

P014

P020

P025

